
EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
 Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch
 EDENTA Etablissement, Industriestrasse 13, FL-9486 Schaanwald – Liechtenstein
 Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

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
Product information and instructions for use for StripFIX:

Product description:

The finishing system consists of four different strip holders with single-sided, perforated diamond strips in three different grit sizes. The system also includes a serrated strip without diamond coating and an ergonomic handle with special grip recess to ensure reliable use. The click-plug system enables easy, quick exchange of the strip holders.

Area of application:

Used for removing cement after cementation of crowns, bridges, onlays or veneers. Used for finishing and polishing interproximal fillings and for removing excess filling material. For dental hygienists to enable removal of discoloration that cannot be accessed with their instruments. In interproximal reposition (IPR) used for reducing interproximal enamel during orthodontic procedures.

Application:

Hold the handle with the company name facing up (1) and attach the strip holder, with the strip facing down, to the tip of the handle (2) and lock into position (3). The strip holder has a mark on the diamond-coated side (6). The swivel-mounted holders can each be locked in 4 positions (4). This ensures that the optimum working position (angle) of the diamond strip is fixed to enable easy access to the anterior and posterior teeth in the oral cavity. Check that the holder is fitted securely before each use. The strip holder can be used with or without the handle. If the strip is used without the handle, it must be secured with a safety chain to protect against aspiration. Hold the strip between the thumb and index or middle finger (5).

Diamond strips:

Insert the strip interdentally using a sawing movement. Face the diamond-coated side towards the filling to be prepared. Move the diamond-coated section over the filling in a grinding movement to remove the excess and adapt the filling anatomically. The strips should be used in descending sequence from medium to super fine grit size.

The colour coding indicates grit size:

SF-M, grit size, medium, blue – 45 µm / 0.13 mm

SF-F, grit size, fine, red – 30 µm / 0.10 mm

SF-SF, grit size, super fine, yellow – 15 µm / 0.08 mm

SF-SS serrated edge without diamond coating - white – 0.04 mm

Important information:

When using a grinding movement, ensure that the edge of the strip does not injure the interdental papillae. After the strip has been inserted, it should be kept closely on the tooth to avoid injuring the papillae. It is essential to avoid bending the strip, as this can result in damage and splitting off of the diamond coating.

Serrated edge:

The serrated edge (7) is used for interdental separation. Once the interproximal has been separated by the serrated edge, replace the strip with a diamond-coated strip and proceed as described above. This avoids injuring the papillae.

Disinfection, cleaning and sterilisation:

The handle and strip holders must be disinfected, cleaned and sterilised before first use on patients and after each use.

Detailed preparation recommendations according to DIN EN ISO 17664  www.edenta.com

Disinfection:

Avoid bending the strips during disinfection, cleaning and sterilisation. Only DGHM (German Society for Hygiene and Microbiology) approved solutions for rotary instruments may be used as a disinfectant solution. The manufacturer's instructions for use of the disinfectant solution should be followed exactly. After the prescribed reaction time, remove the handle and strip holder from the disinfectant solution, rinse under running water and air dry.

Cleaning:

If cleaning is performed in an ultrasonic cleaner, ensure that the strips are not in the same container as other instruments, particularly diamond instruments (possible damage of the strips). Strips can also be cleaned with a brush.

Sterilisation:

Sterilisation can be performed in an autoclave in the usual way. Strips should be sterilised at 134°C; at 2.1 bar; for 5 minutes.

Instrument control:


The strips (strip holders) must be checked for damage before storing. The strip should be discarded if the diamond coating is worn or damaged or the strip is bent.

Storage:

Store dry and protected against recontamination. Protect against bending and other damage.

EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
 Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch
 EDENTA Etablissement, Industriestrasse 13, FL-9486 Schaanwald – Liechtenstein
 Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

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